

KONQUEST **KDT-1201**




QUICK USER GUIDE

KONQUEST KDT-1201

ORAL - ARMPIT - RECTAL THERMOMETER

PRODUCT STRUCTURE



BATTERY REPLACEMENT




Figure 1




Figure 2





Figure 3

1. Replace battery when “” appears in the lower right corner of LCD display.
2. Pull battery cover off as shown in Figure 1.
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2"). Figure 2
4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41, UCC392, or equivalent.
5. Be sure battery is installed with “+” polarity facing up. Figure 3
5. Slide battery chamber back into place and attach cover.

1 DIRECTIONS

[1] Press the On/Off Button next to LCD display. A tone will sound as the screen shows **188.8^{°E}_M**, followed by last recored temperature. After showing the self-test temperature, the thermometer is now in the testing mode.

*Note: If the measured temperature is less than 90°F or 32°C, the LCD will display Lo, if the measured temperature is higher than 109.9°F or 42.9°C, the LCD will display Hi.

[2] Position thermometer in desired location (mouth, rectum, or armpit.)

a) Oral Use: Place thermometer under tongue as indicated by "✓" position shown in Figure 4. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Normal temperature between 96.3°F and 99.1°F (35.7 °C and 37.3 °C)




Figure 4

b) Rectal Use: Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1 cm (less than 1/2 ") into rectum. Normal temperature between 97.2°F and 99.9°F (36.2 °C and 37.7 °C).

c) Armpit Use: Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Normal temperature between 95.4°F and 98.1°F (35.2 °C and 36.7 °C).

[3] The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signaling tone (beep) must be maintained without exception. The measurement continues even after the buzzer notification. The displayed temperature will not change when thermometer is removed from its testing position.

[4] To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.

[5] Store the thermometer in its protective case.

1 CLEANING AND DESINFECTION

Wipe the thermometer with a soft clean cloth.

For stubborn stains, wipe the thermometer with a cloth that has been dampened with water or a neutral detergent solution and then wring thoroughly. Finish by wiping with a soft dry cloth.

For disinfection, 70% Ethanol or Isopropyl alcohol can be used.

Observe the following to prevent damage to the thermometer.

- Do not use benzene, thinner, gasoline or other strong solvents to clean the thermometer.

- Do not attempt to disinfect the sensing section (tip) of the thermometer by immersing in alcohol or in hot water (water over 122°F (50°C)).

- Do not use ultrasonic washing to clean the thermometer.

READ INSTRUCTION MANUAL THOROUGHLY BEFORE USE

KONQUEST